Technical Data for Alicat BIOC-Series Mass Flow Controllers

50 slpm of Full Scale through 500 slpm of Full Scale

Tel: 888-290-6060

Standard Specifications (Contact Alicat for available options.)

| SENSOR PERFORMANCE | | | |
|---|--|--------------------------|--|
| Mass Flow Accuracy at calibration conditions ¹ | ± (0.8% of Reading + 0.2% of Full Scale) | | |
| Repeatability (2σ) | ± (0.2% of Reading + 0.02% of Full Scale) | | |
| Steady State Control Range ² | 0.01% - 100% of Full Scale | | |
| Temperature Sensitivity | Mass Flow Zero and Span Shift: 0.02% Full Scale / °C | | |
| Pressure Sensitivity | Mass Flow Zero and Span Shift: ± (0.08% of Reading + 0.02% Full Scale) / atm | | |
| Operating Temperature Range | -10 to 60°C (consult Alicat for expanded range) | | |
| Temperature Accuracy | ± 0.75°C | | |
| Operating Pressure Full Scale | 160 PSIA (consult Alicat for additional options) | | |
| Pressure Accuracy | Above 1 atm: ± 0.5% of Reading | Below 1 atm: ± 0.07 PSIA | |
| Typical Sensor Response Time | 65 - 255 ms (Adjustable) | | |
| Typical Warm-Up Time | <1s | | |

¹ Stated accuracy is after tare under equilibrium conditions. Extreme gas behavior (especially near state boundaries) can introduce additional flow uncertainties. Consult Alicat if higher accuracy is required.

² Achievable steady state control may be limited by user-configurable PID tuning and process conditions. Dynamic control performance is also limited by control response time, which may vary with the flow rate.

| MECHANICAL | | | |
|---|---|--|--|
| Minimum Operating Pressure | 11.5 PSIA common mode pressure (consult Alicat for lower operating pressures) Differential pressure must exceed model pressure drop, see below for details | | |
| Maximum Operating Pressure | Damage possible above 175 PSIA common mode pressure Damage possible above 75 PSID differential pressure | | |
| Ingress Protection | IP40 (IP66 Option Available) | | |
| Humidity Range | 0 to 95% non-condensing | | |
| Dimensions, pressure drop, weight, and process connection specifications are listed on mechanical drawing pages | | | |

| CONTROL AND COMMUNICATIONS | | | |
|---------------------------------|--|--|--|
| Analog I/O | 0-5 VDC (Serial and Modbus RTU only) | | |
| Digital I/O Options | DeviceNet, EtherCAT, EtherNet/IP, Modbus RTU over RS-232, Modbus RTU over RS-485, Modbus TCP/IP, Profibus, RS-232 Serial, RS-485 Serial | | |
| Electrical Connection | 8 pin M12 or Protocol Dependent | | |
| Power Requirements ³ | 24 VDC, 1,000 mA max | | |
| Data Update Rate ³ | Serial: 40 Hz at 19200 baud Analog: 1000 Hz | | |
| Display Update Rate | 10 Hz | | |
| Analog Signal Accuracy | ± 0.1% of Full Scale additional uncertainty | | |
| Typical Control Response Time | 150 ms to 63% of step change (T63) | | |
| Valve Function | Normally Closed | | |

³ Consult the individual operating bulletins for specific industrial protocol power requirements and data transmission specifications.

| FEATURES | | |
|--|--|--|
| STP Reference Conditions | 25°C and 1 atm (Default), user configurable | |
| NTP Reference Conditions | 0°C and 1 atm (Default), user configurable | |
| Color TFT Display with integrated touchpad | Simultaneously displays Mass Flow, Volumetric Flow, Pressure and Temperature | |
| Gas Select™ | 98 user selectable gases stored internally. Each gas optimized to match NIST's REFPROP 10 gas property calculations across the operating temperature and pressure ranges for highest accuracy. | |
| COMPOSER™ | Allows 20 user definable gas mixes. Up to 5 constituent gases per mix, down to percentages of 0.01% | |

Wetted Materials

| FLOW BODY WETTED MATERIALS | OPTION | VALVE WETTED MATERIALS |
|--|--------|--|
| 316L Stainless Steel, USP VI FDA Certified Viton Elastomers | A | FFKM, 316L Stainless Steel, Elgiloy Super Alloy, Sandvik Super Alloy |
| Each controller has 3 parts: Flow body · Sensor · Valve | В | 302/303/430FR Stainless Steel, Brass, Viton |
| | С | 302/303/304/410 Stainless Steel, Viton, Delrin |
| | OPTION | SENSOR WETTED MATERIALS |
| ASME BPE-2016 Compliance Requires both Valve A and Sensor A | A | 316L Stainless Steel |
| | В | Polyamide, Alumina, Ceramic, Glass, Gold, Silicon, Nylon, Delrin, Heat Cured Epoxy, RTV, Silicone |

DOC-SPECS-BIOC-HIGH 1 REV 0, September 13, 2019


Product Mechanical Drawings and Dimensions



BIOC-Series

0 – 50 slpm

0 - 100 slpm



Flow Range Specific Specifications

| FULL SCALE FLOW MASS CONTROLLER | PRESSURE DROP5 AT FS FLOW (PSID) VENTING TO ATMOSPHERE | APPROXIMATE WEIGHT | MECHANICAL DIMENSIONS ⁶ | PROCESS CONNECTIONS ⁷ |
|------------------------------------|--|-----------------------|---------------------------------------|-------------------------------------|
| 50 slpm | 2.0 | 9.0 lb | 6.4"H x 7.7"W x 2.3"D | %6-18 SAE6 |
| 100 slpm | 3.2 | | | Female |

- 5 Lower Pressure Drops Available, please see our WHISPER-Series mass flow controllers at www.alicat.com/whisper.
- 6 See drawings for metric equvalents.
- 7 Additional process connections available on request. Consult Alicat for more information.


Product Mechanical Drawings and Dimensions





BIOC-Series

0 – 200 slpm

0 - 500 slpm







Flow Range Specific Specifications

| FULL SCALE FLOW MASS CONTROLLER | PRESSURE DROP5 AT FS FLOW (PSID) VENTING TO ATMOSPHERE | APPROXIMATE WEIGHT | MECHANICAL DIMENSIONS ⁶ | PROCESS CONNECTIONS ⁷ |
|------------------------------------|--|-----------------------|---------------------------------------|-------------------------------------|
| 250 slpm | 2.4 | 9.0 lb | 6.7"H x 8.3"W x 2.3"D | 7/8-14 SAE10 |
| 500 slpm | 6.5 | | | Female |

- 5 Lower Pressure Drops Available, please see our WHISPER-Series mass flow controllers at www.alicat.com/whisper.
- 6 See drawings for metric equvalents.
- $7 \quad \text{Additional process connections available on request. Consult Alicat for more information}.$



